

Periodic Acid Schiff (PAS) Stain Kit (Modified Lillie's)

Description: The Periodic Acid Schiff (PAS) Stain Kit is intended for use in histological demonstration of lymphocytes and mucopolysaccharides. The staining pattern of the lymphocytes are helpful in making therapeutic decisions in established cases of lymphocytic leukemia. The PAS reaction in tissue sections is useful for the demonstration of mucopolysaccharides. PAS staining may also be used for the demonstration of fungal organisms in tissue sections.

PAS Positive Material:	Magenta
Nuclei:	Black/Blue


Uses/Limitations: For In-Vitro Diagnostic use only.
 Histological applications.
 Do not use past expiration date.
 Use caution when handling these reagents.

Control Tissue: Kidney
 Intestine
 Liver


Availability/Contents:

<u>Item #</u>	<u>Kit Contents</u>	<u>Storage Conditions</u>
PAQ250	Periodic Acid Solution (250 ml)	2-8° Centigrade
SRF250	Schiff's Solution (250 ml)	2-8° Centigrade
HMM125	Hematoxylin, Mayer's (125 ml)	Room Temperature
BRT125	Bluing Reagent (125ml)	Room Temperature

Precautions: Keep away from open flame.
 Avoid contact with skin and eyes.
 Harmful if swallowed.
 Follow all Federal, State, and local regulations regarding disposal.
 Use in chemical fume hood whenever possible.
 Wear protective clothing.

Storage: 2° C  25° C

**Mixed Storage Conditions.
 Separate Contents.**

 ScyTek Laboratories, Inc.
 205 South 600 West
 Logan, UT 84321
 435-755-9848
 U.S.A.

CE Authorized Representative in Europe
 (Regulatory affairs only)
 EmergoEurope (31)(0) 70 345-8570
 Molsnstraat 15
 2513 BH Hague, The Netherlands

Instructions For Use PAS-1-IFU

Rev. Date: Dec. 20, 2007

Revision: 1

Page 2 of 2


P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

Procedure (Standard):


1. Deparaffinize sections if necessary and hydrate to distilled water.
2. If sections are Zenker-fixed, remove mercuric chloride crystals using iodine and clear with sodium thiosulfate. Rinse in running tap water.
3. Immerse slide in Periodic Acid Solution (PAQ250) for 5 minutes (10 minutes for Kidney, skin and diastase digested liver sections).
4. Rinse slide in 4 changes of distilled water.
5. Immerse slide in Schiff's Solution (SRF250) for 15 minutes (30 minutes for Kidney, skin and diastase digested liver sections).
6. Rinse slide in hot running tap water.
7. Rinse slide in distilled water.
8. Stain slide in Hematoxylin, Mayer's (HMM125) for 2-3 minutes.
9. Rinse slide in running tap water for 2-3 minutes.
10. Apply Bluing Reagent (BRT125) for 30 seconds.
11. Rinse in distilled water.
12. Dehydrate through graded alcohols.
13. Clear, and mount in synthetic resin.

References:

1. Culling CFA, Allison RT, Barr WT.: Cellular Pathology Technique, 4th Edition. Butterworths, Pages 216-220, 1985.
2. Sheenan, D.C., Hrapchak, B.B. Theory and Practice of Histotechnology, 2nd Edition. CV Mosby, Columbus, OH. Pages 164-167, 1980.

Storage: 2° C  25° C

**Mixed Storage Conditions.
Separate Contents.**

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
435-755-9848
U.S.A.

 Authorized Representative in Europe
 (Regulatory affairs only)
EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands