


Twort's Counterstain Kit

Description: *The Twort's Counterstain Kit is designed for histological procedures to provide contrast for easier visualization of tissue sections.*

*Nuclei: Red
Cytoplasm: Light Green*

Uses/Limitations: *Not to be taken internally.
For In-Vitro Diagnostic use only.
Histological applications.
Do not use if reagents become cloudy.
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile.*



Control Tissue: *Any well-fixed tissue.
Cytology Smear.*

Ordering Information and Current Pricing at www.scytek.com


Availability/Contents:


<u>Item #</u>	<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
TCA125	Twort's Neutral Red (Part A)	125 ml	18-25°C
TCB030	Twort's Fast Green (Part B)	30 ml	18-25°C
AAF250	Acetic Alcohol Solution (2%)	250 ml	18-25°C

Precautions: *Avoid contact with skin and eyes.
Flammable.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.*


Prepare the Following Solution Immediately Before Use.

Working Twort's Counterstain Solution: 30 ml Distilled Water
(Add components in the order listed) 9 ml Twort's Neutral Red (Part A)
 2 ml Twort's Fast Green (Part B)

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.


 EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands


Procedure:

1. *Deparaffinize sections if necessary and hydrate to distilled water.*
2. *Stain slide in Working Twort's Counterstain for 5-10 minutes. Use once and discard.*
3. *Rinse quickly in distilled water.*

Note: for greatest retention of stain, air dry, clear and mount at this point. Alcohol and decolorizer solutions can quickly remove stains.


4. *Dip once in absolute alcohol.*
5. *Decolorize quickly and carefully in Acetic Alcohol Solution.*
6. *Dip two times in absolute alcohol.*
7. *Clear, and mount in synthetic resin.*

Storage: 18° C  25° C



ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.

CE 

 EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands