

Alkaline Phosphatase Enhancer

Description: A quick and effective solution for enhancing the signal of chromogens that utilize Alkaline Phosphatase as the enzyme. ScyTek's Alk-Phos kits do not need this solution for maximum signal but other systems may benefit from its use.

Uses/Limitations: Not to be taken internally. Histological applications. Do not use if reagents become cloudy. Do not use past expiration date. Use caution when handling reagents. Non-Sterile.

Storage: Room temperature (18-25°C)

Availability:	Item #	Volume
	APE125	125 ml
	APE500	500 ml
	APE999	1000 ml

Precautions: Avoid contact with clothes and exposed skin. May be harmful if swallowed. Follow all federal, state, and local regulations regarding disposal.

Procedure:

1. Once tissue sections have been incubated with alkaline phosphatase, wash them with buffer thoroughly.
2. Wipe the glass to remove excess buffer.
3. Incubate for 1 minute at room temperature.
4. Drain excess alkaline phosphatase enhancer and add substrate/chromogen without washing.

Storage: 18° C



25° C



ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.



EC REP

Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands