

Buffered Glycerol Mounting Media

Description: Buffered Glycerol Mounting Media is a non-drying formulation that is useful in a variety of procedures, such as immunofluorescence, immunohistochemistry, special stains, etc. Coverslips may be easily removed by soaking in water.
No heating is required before use.
Contains no Anti-Fade components.
pH 8.0

Availability/Contents:	Item #	Volume
	GMM003	3 ml Dropper
	GMM030	30 ml Dropper
	GMM500	500 ml
	GMM999	1000 ml

Uses/Limitations: Not to be taken internally.
For In-Vitro Diagnostic use.
Histological applications.
Do not use if reagents become cloudy.
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile.



Ordering Information and Current Pricing at www.scytek.com


Precautions: Avoid contact with skin and eyes.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.

Procedure:

1. Apply 1-2 drops of Buffered Glycerol Mounting Media to microscope slide and coverslip.
2. If desired, seal edges with nail polish or synthetic resin mounting media.
3. For immunofluorescent procedures, store slides at 2-8°C in the dark.

Note: Coverslips may be removed by soaking in water.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.

CE 

EC REP

Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands