


## Gram Stain Kit

**Description:** The Gram Stain Kit is intended for the demonstration and differentiation of Gram-positive and Gram-negative bacteria. The stain is named after the Danish scientist Hans Christian Gram.

Gram Positive Bacteria: Blue  
 Gram Negative Bacteria: Red  
 Other Tissue: Yellow  
 Nuclei: Red

**Uses/Limitations:** Not to be taken internally.  
 For In-Vitro Diagnostic use only.  
 Histological applications.  
 Do not use if reagents become cloudy.  
 Do not use past expiration date.  
 Use caution when handling reagents.  
 Non-Sterile.



**Control Tissue:** Any well fixed tissue section.  
 Air dried smear.

**Ordering information regarding individual components on back page!**


**Kit Contents:**


<u>Item #</u>	<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
GVS125	Gentian Violet Solution	125 ml	18-25 °C
LIS125	Lugol's Iodine Solution	125 ml	18-25 °C
GDS125	Gram's Decolorizer Solution	125 ml	18-25 °C
CFX125	Carbol Fuchsin Counterstain	125 ml	18-25 °C
TZQ125	Tartrazine Solution	125 ml	18-25 °C

**Precautions:** Avoid contact with skin and eyes.  
 Harmful if swallowed.  
 Follow all Federal, State, and local regulations regarding disposal.


**Procedure:**

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Apply adequate Gentian Violet Solution to completely cover tissue section and incubate for 1 minute.
3. Rinse slide in distilled water to remove excess stain.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.  
 205 South 600 West  
 Logan, UT 84321  
 U.S.A.

 EmergoEurope (31)(0) 70 345-8570  
 Molsnstraat 15  
 2513 BH Hague, The Netherlands

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 – Tel. (435) 755-9848 - Fax (435) 755-0015 - [www.scytek.com](http://www.scytek.com)


4. Apply adequate Lugol's Iodine Solution to completely cover tissue section and incubate for 1 minute.
5. Rinse slide in gently running tap water to remove excess iodine.
6. Place slide in Gram's Decolorizer until color no longer bleeds off section.
7. Rinse slide quickly in gently running tap water.
8. Apply adequate Carbol Fuchsin to completely cover tissue section and incubate for 1-2 minutes.
9. Rinse slide quickly in gently running tap water to remove excess stain.
10. Apply Tartrazine Solution and incubate for 15 seconds.
11. Rinse slide 1 time in absolute alcohol.
12. Dehydrate slide quickly in 3 changes of absolute alcohol.
13. Clear in 2 changes of xylene or xylene substitute, and mount in synthetic resin.

#### References:

1. Sheehan, DC., Hrapchak, BB. Theory and Practice of Histotechnology; 1980, page 235.
2. Su, R.J., Wang, P. Role of Gram stain in microbiological laboratories with limited resources. Reviews in Medical Microbiology. July 2011, Volume 22, Issue 3: pages 41-44. Doi: 10.1097/MRM.0b013e3283478a08.
3. Marira, J., Surekha, Y, Asangi, K.S., Suresh, B.S., Ramesh, S. Sputum Gram Stain Assessment in Relation to Sputum Culture for Respiratory Tract Infections in a Tertiary Care Hospital. Journal of Clinical and Diagnostic Research. December 2011, Volume 5(8): pages 1699-1700.


#### Bulk Reagent Ordering Information and Current Pricing at [www.scytek.com](http://www.scytek.com)

Description:	Catalog #	Volume
Gentian Violet Solution	GVS125	125 ml
	GVS500	500 ml
	GVS999	1000 ml
Lugol's Iodine Solution	LIS125	125 ml
	LIS500	500 ml
	LIS999	1000 ml
Gram's Decolorizer	GDS125	125 ml
	GDS500	500 ml
	GDS999	1000ml
Carbol Fuchsin Counterstain	CFX125	125 ml
	CFX500	500 ml
	CFX999	1000 ml
Tartrazine Solution	TZQ125	125 ml
	TZQ500	500 ml
	TZQ999	1000 ml

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
U.S.A.



 EmergoEurope (31)(0) 70 345-8570  
Molsnstraat 15  
2513 BH Hague, The Netherlands