

Picro-Sirius Red Stain Kit (For Cardiac Muscle)

Description:

The Picro-Sirius Red Stain Kit (For Cardiac Muscle) is intended for use in the histological visualization of thin septa and collagen fibers. This modification of the Picro-Sirius Red stain eliminates the yellow cytoplasmic staining that can obscure thin collagenous septa. Using this procedure easily allows viewing of collagenous septa as thin as 0.2-0.5 microns. The PSR stain may be viewed using standard light microscopy or polarized light resulting in birefringence of the collagen fibers.

Light Microscopy

Collagen: Red
 Septa: Red
 Cytoplasm: Colorless to Slightly Yellow

Polarized Light Microscopy


Collagen: Yellow-Orange and Green Birefringence

Uses/Limitations:

Not to be taken internally.
 For In-Vitro Diagnostic use only.
 Histological applications.
 Do not use if reagents become cloudy.
 Do not use past expiration date.
 Use caution when handling reagents.
 Non-Sterile.

Control Tissue:

Cardiac Muscle
 Uterus
 Lung
 Kidney



Ordering information regarding individual components on back page!

Availability/Contents:

<u>Item #</u>	<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
PQA250	Phosphomolybdic Acid Solution (0.2%)	250 ml	18-25°C
SRS250	Picro-Sirius Red Solution	250 ml	18-25°C
AAD250	Acetic Acid Solution (0.5%)	250 ml	18-25°C

Precautions:

Avoid contact with skin and eyes.
 Harmful if swallowed.
 Follow all Federal, State, and local regulations regarding disposal.
 Use in chemical fume hood whenever possible.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
 205 South 600 West
 Logan, UT 84321
 U.S.A.

CE 

EC REP
 Emergo Europe
 Prinsessegracht 20
 2514 AP The Hague, The Netherlands

Procedure:

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Apply adequate Phosphomolybdic Acid Solution (0.2%) to completely cover tissue section and incubate for 1-5 minutes.
3. Dip slide one time in distilled water.
4. Apply adequate Picro-Sirius Red Solution to completely cover tissue section and incubate for 60-90 minutes.
5. Rinse slide quickly in two changes of Acetic Acid Solution (0.5%).
6. Rinse slide using absolute alcohol.
7. Dehydrate in 2 changes of absolute alcohol, clear, and mount in synthetic resin.

References:

1. Puchtler H., Waldrop F.S., Valentine L.S. Polarization microscopic studies of connective tissue stained with picro-sirius red FBA. Beitr Path. 1973; 150, pages 174-187.
2. Junqueira L.C.U., Bignolas G., Brentani R.R. Picrosirius staining plus polarization microscopy, a specific method for collagen detection in tissue sections. Histochemistry J. 1979, 11, pages 447-455.
3. Dolber, P.C., Spach, M.S. Picrosirius red staining of cardiac muscle following phosphomolybdic acid treatment. Stain Technology, January 1987, Volume 62(1): pages 23-26.
4. Whittaker, P., Kloner, R.A., Boughner, D.R., Pickering, J.G. Quantitative assessment of myocardial collagen with picrosirius red and circularly polarized light. Basic Research in Cardiology. 1994, Volume 89, Number 5, pages 397-410. DOI: 10.1007/BF00788278
5. Whittaker P. Polarized light microscopy in biomedical research. Microscopy and Analysis 1995; 44, pages 15-17.

Bulk Reagent Ordering Information and Current Pricing at www.scytek.com

Description:	Catalog #	Volume
Phosphomolybdic Acid Solution (0.2%)	PQA250	250 ml
	PQA500	500 ml
	PQA999	1000 ml
Picro-Sirius Red Solution	SRS250	250 ml
	SRS500	500 ml
	SRS999	1000 ml
Acetic Acid Solution	AAD250	250 ml
	AAD500	500 ml
	AAD999	1000 ml

Storage: 18° C  25° C



ScyTek Laboratories, Inc.
 205 South 600 West
 Logan, UT 84321
 U.S.A.



EC REP

Emergo Europe
 Prinsessegracht 20
 2514 AP The Hague, The Netherlands