


UltraTek Complete HRP Anti-Polyvalent (DAB) Staining System

Description: The UltraTek Complete HRP Anti-Polyvalent (DAB) Staining System provides unmatched sensitivity with incubation times of 10 minutes each for the Link Antibody and Enzyme Label

Species of Origin:	Goat
Antigen Specificity:	Anti-Polyvalent (Mouse, Rat, Rabbit and Guinea Pig)
Preadsorbed Against:	Human
Enzyme Conjugate:	Peroxidase
Chromogen Substrate:	3,3'-Diaminobenzidine (DAB)

Uses/Limitations: Not to be taken internally.
For Invitro Diagnostic Use.
Histological applications.
Do not use if reagents become cloudy.
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile.




Control Tissue: Any well-fixed tissue section.
Frozen tissue section.
Cytocentrifuge preparation.


Ordering Information and Current Pricing at www.scytek.com

Kit Contents:	<u>Item #</u>	<u>Volume</u>
	Peroxide Block	1 x 15 ml
	Super Block	1 x 15 ml
	UltraTek Anti-Polyvalent	1 x 15 ml
	UltraTek HRP	1 x 15 ml
	DAB Chromogen	1 x 3 ml
	DAB Substrate (High Contrast)	5 x 5 ml vials
	Hematoxylin, Mayer's (Lillie's Mod.)	1 x 15 ml
	Bluing Reagent	1 x 15 ml


Storage: Store at 2-8°C.

Precautions: Avoid contact with skin and eyes.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.

Storage: 2° C  8° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.



 EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands

Procedure:

1. Deparaffinize and rehydrate tissue section.
2. To reduce nonspecific background staining due to endogenous peroxidase, incubate slide in Peroxide Block for 10 minutes.
3. Wash 2 times in buffer.
4. If required, incubate tissue in digestive enzyme.
5. Wash 4 times in buffer.
6. Apply Super Block (blue cap), and incubate for 5 minutes at room temperature to block nonspecific background staining.
Note: Do not exceed 10 minutes or there may be a reduction in desired stain.
7. Wash 1 time in buffer.
8. Apply primary antibody and incubate according to manufacturer's protocol.
9. Wash 4 times in buffer.
10. Apply UltraTek Anti-Polyvalent (yellow cap), and incubate for 10 minutes at room temperature.
11. Wash 4 times in buffer.
12. Apply UltraTek HRP (red cap), and incubate for 10 minutes at room temperature.
13. Rinse 3 times in buffer.

WARNING: DAB is a suspected carcinogen. Handle with care and dispose of according to all regulations.

14. Add 5 drops (200-250ul) DAB Chromogen to one vial of DAB Substrate High Contrast, mix by swirling and apply to tissue for 5 minutes.
15. Rinse 1 time in distilled or deionized water.
16. Apply DAB Chromogen/Substrate mixture and incubate for a second 5 minute period.
17. Rinse 3 times in distilled or deionized Water.
18. Apply Hematoxylin, Mayer's (Lillie's Modification) and incubate for 1-2 minutes.
19. Rinse 3 times in distilled or deionized water.
20. Apply Bluing Reagent for 5-10 seconds.
21. Rinse immediately in distilled or deionized water.
22. Dehydrate slides and clear in xylene or xylene substitute.
23. Coverslip using a permanent mounting media.

Storage: 2° C  8° C ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A. IVDEmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands

-Troubleshooting Guide-

Overstaining:

1. Concentration of the primary antibody was too high or the incubation time was too long.
2. Temperature during incubation was too high.
3. Incubation time with link antibody or streptavidin/enzyme label was too long.

Nonspecific Background Staining:

1. Rinsing between steps was inadequate.
2. Tissue was allowed to dry with reagents on.
3. Folds in tissue trapped reagents.
4. Tissue contains endogenous peroxidase.
5. Tissue contains endogenous biotin.
6. Antigen migrated in tissue.
7. Excessive tissue adhesive on slides.
8. Inadequate blocking with protein block.


Weak Staining:

1. Primary antibody concentration was too low or incubation time was too short.
2. Reagents are past their expiration date.
3. Inadequate removal of wash water between steps, resulting in dilution of reagents.
4. Counterstain or mounting media were incompatible and dissolved the chromogen reaction product.
5. Room temperature was excessively cool.
6. The primary antibody does not recognize an antigen that survives fixation and embedding in high enough amounts.
7. Excessive incubation with protein block (Super Block or normal serum).


No Staining:

1. Steps were inadvertently left out.
2. There is no antigen in the tissue.
3. The primary antibody is not of mouse, rat, rabbit or guinea pig origin.
4. Chromogenic substrate has been replaced with another that is not intended for use with peroxidase.
5. One or more components of the kit have been inactivated.

Storage: 2° C  8° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.

 EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands