

## Normal Antibody Diluent (Phosphate Buffered)

**Description:** Normal Antibody Diluent is prepared for use in immunochemical procedures. Antibodies can be diluted directly in this reagent and can be stored for up to 18 months when refrigerated (extremely high dilutions will result in reduced shelf life). Dilute primary antibodies, or biotinylated antibodies in this reagent.

Appearance: Transparent liquid, slightly yellow. Discard if reagent becomes cloudy.

pH: 7.4±0.1

**Uses:** Immunochemical procedures.  
Antibody diluent.

**Limitations:** Do not use past expiration date.

<b>Availability:</b>	<b>REF</b>	<b>Volume</b>
	ABB125	125 ml
	ABB500	500 ml
	ABB999	1000 ml

**Storage:** Store at 2-8°C. Product is stable for 24 months from the date of manufacture.

**Procedure:**

1. Dilute antibody directly in this reagent.
2. Store diluted antibody at 2-8°.

**Precautions:** Do not pipette by mouth.  
Avoid contact with skin and eyes.  
Observe all federal, state and local environmental regulations regarding disposal.

Storage: 2° C



8° C



ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
435-755-9848  
U.S.A.



Authorized Representative in Europe  
(Regulatory affairs only)  
EmergoEurope (31)(0) 70 345-8570  
Molsnstraat 15  
2513 BH Hague, The Netherlands