

Acid Fuchsin (0.5% Aqueous)

Description: This stain is for histology applications.

Results:	Nuclei	Red
	Fibroglia and Myoglia	Red
	Neuroglia Fibrils	Red
	Axis Cylinders and Fibrin	Red
	Nucleoli	Red

Availability:

<u>Item #</u>	<u>Volume</u>
AF1500	500 ml
AF1999	1 Liter

Procedure:

1. Remove paraffin and re-hydrate as usual.
2. Stain 1 - 5 minutes in Acid Fuchsin solution.
3. Transfer directly to 95% alcohol. Change alcohol several times to remove excess stain.
4. Dehydrate in absolute alcohol, clear and mount.