

# Polyclonal Antibody to Human Lambda Light Chain

## Description:

- Immunogen:** Lambda light chain isolated from pooled urine of patients with Bence Jones proteinuria.
- Format:** This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with free as well as bound lambda light chains. Contaminating antibodies have been removed by solid phase absorption.

## Uses/Limitations:

Immunohistochemistry  
For Research Use Only.  
Do not use past expiration date.

## Storage:

2-8° Centigrade.

## Procedure:

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

## Precautions:

Contains Sodium Azide as a preservative.  
Do not pipette by mouth.