

## Polyclonal Antibody to Desmin

### Description:

Immunogen:	Desmin isolated from chicken gizzard.
Host:	Rabbit
Format:	Whole rabbit antiserum against chicken smooth muscle desmincontaining sodium azide as preservative.
Specificity:	This antibody reacts with both smooth and striated muscle cells. This antibody can be used for typing of soft tissue sarcomas.
Titer/Working Dilution:	This antibody may be diluted to a titer of 1:100-1:250 in an immunohistochemistry application. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.

**Uses/Limitations:** Immunohistochemistry  
For Research Use Only.  
Do not use past expiration date.

**Storage:** 2-8° Centigrade.

**Procedure:** We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user.

**Precautions:** Contains Sodium Azide as a preservative.  
Do not pipette by mouth.