

Zinc Formalin

Description: Zinc Formalin is a tissue fixative for routine procedures.

Contents: Water, Formalin (10%), Zinc Chloride, Glacial Acetic Acid.

pH: 4.0 (Adjust if desired).

Storage: Store at room temperature. Do not use beyond the expiration date.

Procedure:

1. Trim specimen to smallest practical size.
2. Cover tissue with an amount of fluid equal to 20 times the specimen volume.
3. Fixation times range from 5 to 48 hours depending on specimen size.

Precautions: Poison. Do Not Take Internally.