

Polyclonal Antibody to Human Prostate Specific Antigen (PSA)

Description:

- Immunogen:** Prostate specific antigen isolated from human seminal plasma.
- Format:** This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with PSA. The staining was observed in the prostatic ductal epithelium of normal and neoplastic prostate. No staining with any other cell type has been observed.

Uses/Limitations: Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage: 2-8° Centigrade.

Procedure: We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user.

Precautions: Contains Sodium Azide as a preservative.
Do not pipette by mouth.