

Specification Sheet

Monoclonal Mouse Antibody to Human Prostatic Acid Phosphatase (PSAP)

Catalog No.:	A00041
Immunogen:	Prostatic acid phosphatase purified from human seminal plasma.
Clone:	PASE/4LJ.
Isotype:	IgG1.
Format:	This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
Staining Protocol:	We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.
Specificity:	This antibody is specific against PSAP. This antibody can be used to identify the prostatic origin of metastatic tumor, for example metastatic carcinoma in bone, rectum and lymph node.
Storage:	Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.
References:	i) Haines et al. Biochem Soc Trans 15: 1179, 1987. ii) Haines et al. British J Cancer 60: 887, 1989.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

ScyTek Laboratories, Inc. will not be held responsible for patent infringement or other violation that may occur with the use of our product.