

Polyclonal Antibody to Human IgM

Description:

- Immunogen:** IgM isolated from a pool of normal human sera.
- Format:** This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with the mu-chains of human IgM. Nonspecific antibodies have been removed by solid-phase absorption.

Uses/Limitations:

Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage:

2-8° Centigrade.

Procedure:

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

Precautions:

Contains Sodium Azide as a preservative.
Do not pipette by mouth.