

Polyclonal Antibody to Human Myoglobin

Description:

- Immunogen:** Myoglobin isolated from human heart.
- Format:** This antibody has been pre-titrated and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with human myoglobin. The antibody stains strongly with the skeletal and cardiac muscles. No non-specific staining with other tissues has been observed.

Uses/Limitations:

Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage:

2-8° Centigrade.

Procedure:

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

Precautions:

Contains Sodium Azide as a preservative.
Do not pipette by mouth.