

## Enzyme Label Diluent

**Description:** Enzyme Label Diluent is prepared for use in immunochemical procedures. Enzyme labeled avidin/streptavidin can be diluted directly in this reagent and can be stored for 18 months (immunohistochemical dilutions) when refrigerated. High dilutions result in decreased shelf life. This product is subjected to 0.2 micron filtration.

Appearance: Transparent liquid, slightly yellow. Discard if reagent becomes cloudy.

pH: 7.6±0.1

**Availability/Contents:**

<u>Item #</u>	<u>Volume</u>
ABI125	125 ml
ABI500	500 ml
ABI999	1000 ml

**Storage:** Store at 2-8°C. Product is stable for 24 months from the date of manufacture.

**Procedure:** Dilute enzyme labeled streptavidin/avidin/antibody directly into this product.

**Precautions:** Do not pipette by mouth. Avoid contact with skin and eyes. Observe all federal, state and local environmental regulations regarding disposal.