

Polyclonal Antibody to Human IgG

Description:

- Immunogen:** IgG isolated from a pool of normal human sera.
- Format:** This antibody has been pretitrated and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with gamma-chains of human IgA. Nonspecific antibodies have been removed by solid-phase absorption.

Uses/Limitations:

Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage:

2-8° Centigrade.

Procedure:

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Enzymatic treatment with proteolytic enzyme may be required with certain tissues.

Precautions:

Contains Sodium Azide as a preservative.
Do not pipette by mouth.