

Polyclonal Antibody to Human Calcitonin

Description:

Immunogen:	Synthetic.
Host:	Rabbit
Format:	This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
Specificity:	This antibody reacts with the mu-chains of human IgM. Nonspecific antibodies have been removed by solid-phase absorption

Uses/Limitations: Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage: 2-8° Centigrade.

Procedure: We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

Precautions: Contains Sodium Azide as a preservative.
Do not pipette by mouth.