

Normal Antibody Diluent (TRIS Buffered)

Description: Normal Antibody Diluent (TRIS Buffered) is prepared using the highest quality components for use in immunochemical procedures. Antibodies can be diluted directly in this reagent and stored for up to 18 months when refrigerated (extremely high dilutions will result in reduced shelf life).

Availability:

<u>Item #</u>	<u>Volume</u>
ADT125	125 ml
ADT500	500 ml
ADT999	1000 ml

Storage: Store at 2-8°C.

Procedure: Dilute antibody directly in this reagent.

Precautions: Do not pipet by mouth. Avoid contact with skin and eyes. Observe all federal, state and local environmental regulations regarding disposal.