

## Safranin O Solution


**Description:** Safranin O is commonly used for counterstaining nuclei red.

**Availability:**

<u>Item #</u>	<u>Volume</u>
SOH500	500 ml
SOH999	1000 ml

**Procedure:**

1. Re-hydrate tissue as usual.
2. Stain tissue section for 5 minutes with Safranin O Solution.
3. Rinse in 95% alcohol 1 minute followed by 1 minute in absolute alcohol.
4. Clear in xylene, and mount in permanent mounting medium.

Storage: 18° C  25° C



ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
435-755-9848  
U.S.A.



Authorized Representative in Europe



(Regulatory affairs only)  
EmergoEurope (31)(0) 70 345-8570  
Molsnstraat 15  
2513 BH Hague, The Netherlands