

## Tween® 20 (Polysorbate 20)

**Description:** Tween® 20 (Polysorbate 20) is used in a broad range of biotechnical applications.

**Uses/Limitations:** Not to be taken internally.  
For Research Use Only.  
Histological applications.  
Do not use if reagent becomes cloudy.  
Do not use past expiration date.  
Use caution when handling reagents.  
Non-Sterile.

**Storage:** Store at 18-25°C.


**Precautions:** Avoid contact with skin and eyes.  
Harmful if swallowed.  
Follow all Federal, State, and local regulations regarding disposal.

**Procedures:**

1. When used with Phosphate Buffered Saline or Tris Buffered Saline as a washing agent, add Tween® 20 at a concentration of 0.05% v/v.
2. When used as a blocking agent, add Tween® 20 at a concentration 0.05 – 0.5% v/v.

Tween® 20 is a registered trademark of Croda Americas.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
U.S.A.

CE

EC REP

Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague, The Netherlands